

Micro DSAEK. Macro Benefits.

New Preloaded
LEITR DSAEK 3.0

Priya M. Mathews, M.D., M.P.H.

TransplantREADY™

Prepared and preloaded for
small incision, no-touch delivery.

*“Easiest, fastest DSAEK I’ve performed since
2007; the graft basically opened by itself.”*

—Thomas A. Planchard M.D.

*“LEITR’s DSAEK 3.0 is simple and straightforward;
the technique is elegant—I can’t see performing
DSAEK any other way.”*

—Aarup A. Kubal, M.D.



1410 N 21st Street, Tampa, FL 33605-5313
813.289.1200 | lionseyeinstitute.org



Register now to evaluate
LEITR DSAEK 3.0.

Scan QR Code or visit:
LionsEyeInstitute.org/dsaek-3/

Patent pending. DSAEK graft stained for illustration purposes.

Access the Benefits.

Optimal Convenience,
Cosmesis & Durability.

OptiGraft™

Sterile Ophthalmic Allografts

- Efficient process designed to ensure safety, quality and affordability
- Ready-to-use tissue without re-hydration or rinsing
- Stores at room temperature
- Easy-access container design

"We used two OptiGraft half-moon patch grafts today. Loved them! Packaging is so much better than that of the competitor."
—Aarup Kubal, M.D.

*David D. Moore, CEBT Laboratory
Operations Manager at LEITR*



1410 N 21st Street, Tampa, FL 33605-5313
813.289.1200 | lionseyeinstitute.org



Explore the full line of sterile graft types and sizes.

Scan QR Code or visit:
LionsEyeInstitute.org/OptiGraft/